



## Product Sheet ~ Elementar CHNRapid

### ~ Aluminium boats

Cat #	Description	OEM #
D5132	Aluminium Weighing Boats 12 x 4 x 4mm pack of 1000	-
D5133	Aluminium Weighing Boats 12 x 6 x 6mm pack of 500	-

### ~ Aluminium capsules for liquids

Cat #	Description	OEM #
D4111	Aluminium Capsules Smooth Wall 7 x 3.2mm pack of 250	05 000 814
D4110	Aluminium Capsules Smooth Wall 9 x 3.5mm pack of 250	05 000 815

### ~ Silver boats

Cat #	Description	OEM #
D5026	Silver Weighing Boats 12 x 6 x 6mm pack of 100	22 133 212

### ~ Silver capsules for liquids

Cat #	Description	OEM #
D4102	Silver Capsules Smooth Wall 5.5 x 3.5mm pack of 250	05 000 816
D4101	Silver Capsules Smooth Wall 9 x 3.5mm pack of 100	-

### ~ Tin boats

Cat #	Description	OEM #
D5022	Tin Weighing Boats 12 x 4 x 4mm pack of 100	-
D5130	Tin Weighing Boats 12 x 4 x 4mm pack of 1000	22 137 418
D5023	Tin Weighing Boats 12 x 6 x 6mm pack of 100	-
D5131	Tin Weighing Boats 12 x 6 x 6mm pack of 500	22 137 419
D5140	Tin Weighing Boats Heavy Weight 12 x 4 x 4mm pack of 1000	-

### ~ Tin capsules for liquids

Cat #	Description	OEM #
D4121	Smooth Wall Tin Capsules Flat Base 11 x 4.75mm pack of 250	05 000 481
D4033	Smooth Wall Tin Capsules Flat Base 5.5 x 3.5mm pack of 100	-
D4001	Smooth Wall Tin Capsules Flat Base 6 x 3mm pack of 100	-
D4131	Smooth Wall Tin Capsules Flat Base 6 x 3mm pack of 250	03 951 620
D4035	Smooth Wall Tin Capsules Flat Base 9 x 3.5mm pack of 100	-
D4132	Smooth Wall Tin Capsules Flat Base 9 x 3.5mm pack of 250	03 951 619

### ~ Tin capsules/pans for solids

Cat #	Description	OEM #
D1011	Tin Capsules Pressed Standard Weight 20 x 8mm pack of 100	-

### Chemical Reagents

Cat #	Description	OEM #
B1089	Alumina (Aluminium Oxide) Wool Fine 05000086 20gm	05 000 086
B1219	Cerium Dioxide Granular 03679915 100gm	03 679 915
B1001	Cobaltous/ic Oxide Granular 0.85 to 1.7mm 25gm	05 001 070

B1122	Copper Oxide Wires Fine Wires 4 x 0.5mm 250gm	05 001 039
B1091M	Copper Wires 3x100gm Pack Fine Wires Reduced 4 x 0.5mm 300gm	-
B1015M	Copper Wires Fine Wires Reduced 4 x 0.5mm 100gm	05 000 699
B1083M	Copper Wires Pound Pack Fine Wires Reduced 4 x 0.5mm 454gm	-
B1161	EMADrone Magnesium Perchlorate 501-171 454gm	03 679 901
B1115	EMASorb A Granular 0.8 to 1.6mm 14 to 25 mesh CO2 Absorbant 200gm	03 679 900
B1043	High Purity Carbon Granular 0.3 to 0.85mm 5gm	03 679 910
B1207	Lead(II) Chromate Granular 0.85 to 1.7mm 03679914 250gm	03 679 914
B1029	Lead(II) Chromate Granular 0.85 to 1.7mm 50gm	-
B1204	Manganese Dioxide Sulphur Reduction in U-Tube 05000880 10gm	05 000 880
B1276	Quartz wool 15 g	-
B1197	Sicapent 03960352 (2.8L) 500g	03 960 352
B1214	Silica (Quartz) Chips Granular 1 to 2mm 03679907 100gm	03 679 907
B1215	Silica (Quartz) Chips Granular 2 to 4mm 05000900 50g	05 000 900
B1216	Silica (Quartz) Chips Granular 2 to 4mm 100g	-
B1208	Silica (Quartz) Wool Very Fine 10gm	03 679 908
B1066	Silica Gel 03960157 1kg	03 960 157
B1231	Silver Wool 0.05mm Fine Wire 50gm	22 131 365
B1217	Tungsten (VI) Oxide Granular 11.02-0008 60gm	11.02-0008

### O-rings and operational accessories

Cat #	Description	OEM #
E2505	Gloves 05000094 Pair	05 000 094
E1088	O Ring 03 654 772	03 654 772
E1081	O Ring Nitrile Rubber 03654615	03 654 615
E1085	O Ring Silicone Rubber 03653118	03 653 118
E1086	O Ring Silicone Rubber 03654505	03 654 505
E1016	O Ring Silicone Rubber 03654603 pack of 10	03 654 603
E1066	O Ring Silicone Rubber 03654629	03 654 629
E1082	O Ring Silicone Rubber 03654701	03 654 701
E1055	O Ring Silicone Rubber 03654814	03 654 814
E1090	O Ring Silicone Rubber 03679941	03 679 941
E1062	O Ring Silicone Rubber 05000043	05 000 043
E1083	O Ring Silicone Rubber 05000248	05 000 248
E1073	O Ring Viton Rubber 03654627	03 654 627
E1057	O Ring Viton Rubber 05000369	05 000 369
E1079	Sealing Ring 03667990	03 667 990
E1080	Sealing Ring 03667991	03 667 991
E2039	Snoop Leak Detection Fluid	-
E2240	Vacuum Grease 03968006 100gm	03 968 006

### Pure Chemical Standards

Cat #	Description	OEM #
B2015	4-Nitroaniline OAS 1gm See Certificate 2755 Expiry 19-May-14	-
B2114	Acetanilide OAS See Certificate 151853 Expiry 27 Apr 15 5gm	03 679 942

## Product Sheet ~ Elementar CHNRapid

B2047	Benzoic Acid OAS See Certificate 3003 Expiry 03 Dec 14 5gm	05 000 400
B2147	Sulphanilic Acid OAS 5gm Certificate 152291 Exp 06 may 15	03 001 118

### Quartz, Glassware and Furnace tubes

Cat #	Description	OEM #
C3030	Ash Finger Ceramic Rapid N 05000087	05 000 087
C1020	Cap Kit for Threaded Joints Caps & PTFE Backed Silicone Seals pack of 2	-
C3023	Combustion Tube Insert Transparent Silica CHN Rapid	-
C3004	Combustion Tube Rapid CHN 22156202/3	22 156 202/3
C3008	Drying Guard Tube Borosilicate Elementar 22 160 200/4	22 160 200/4
C3007	Purification U Tube Borosilicate Glass Rapid CHN	-
C3015	Purification U Tube Vario-EL 22160301/4	22 160 301/4
C3020	Pyrolysis Tube Insert Vario-EL 22133211/4	22 133 211/4
C3019	Pyrolysis Tube Vario-EL 22133210/3	22 133 210/3
D6001	Receptacle Insert CHN Rapid 22127313/4	22 127 313/4
C3005	Reduction Tube Rapid CHN 22127245/3	22 127 245/3

**Placing your order:** We are happy to accept your order by email, telephone, fax or post. We welcome all orders, large or small, with no minimum order requirement. Visa and Mastercard payments are accepted or we can open an account for you – please enquire. Most orders are despatched from stock within 24 hours.

**Quality:** All products are manufactured in our ISO9001:2000 approved facility and are covered by our 'Total Satisfaction Guarantee'.

**Can't see what you need?** We offer thousands of products for all elemental analysers. If you can't see what you need here the chances are we will have the product available. Just contact our friendly sales team who will track down what you need.

**OEM Part numbers:** The use of OEM part numbers is intended for convenience only and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and suitable for use in the stated instrument. All trade marks acknowledged.