



high quality consumables for Nitrogen/Protein elemental analysis



Product Sheet ~ Thermo NA1112Protein

Chemical Reagents

Cat #	Description	OEM #
B3017	Chromosorb W 30-60 mesh Acid Washed 33837530 25gm	33837530
B1091	Copper Wires 3x100gm Pack Fine Wires Reduced 4 x 0.5mm 300gm	-
B1016	Copper Wires Coarse Wires Reduced 6 x 0.65mm 100gm	-
B1046	Copper Wires Coarse Wires Reduced 6 x 0.65mm 50gm	33835312
B1199	Copper Wires Pound Pack Coarse Wires Reduced 6 x 0.65mm 454gm	-
B1056	Glass Wool Very Fine 50gm	-
B1101	Magnesium Perchlorate Granular 0.7 to 1.2mm H2O Absorbant 100gm	-
B1047	Molecular Sieve 3A 1/16 pellets 33801801 100gm	33801801
B1103	Oxidation Catalyst NA2000/NA2100 Copper Oxide/Platinised Alumina 40gm	33840000
B1102	Quartz (Silica) Wool Fine 502-177 50gm	-
B1276	Quartz wool 15 g	33822205
B1252	Quartz Wool Discs 22mm Diameter - Pack of 25	-
B1208	Silica (Quartz) Wool Very Fine 10gm	-
B1017	Silica (Quartz) Wool Very Fine 5gm	33822200
B1189	Soda Lime Granular 1 to 1.25mm Self-indicating CO2 Absorbant 200gm	33835230
B1211	Soda Lime Granular 1 to 2.5mm Self-indicating CO2 Absorbant 100gm	33835235
B1092	Soda Lime Granular 1 to 2.5mm Self-indicating CO2 Absorbant 800gm	-

Crucibles and Boats

Cat #	Description	OEM #
D6050	Stainless Steel Crucible Flash EA 1112 25205010	25205010

O-rings and operational accessories

Cat #	Description	OEM #
F6201	25mm bottom fittings complete with interconnecting tubing and fittings FlashEA1112	-
F4506	Aluminium Ferrule 2mm 29034045 pack of 10	29034045
E6051	Bottom O-Ring For Alloy Reactor Flash EA1112 29020649 pack of 5	29020649
E2035	Column Fitting 6 to 2mm Stainless Steel 34715401	34715401
E2016	Forceps Stainless Steel Curved Points	20500500
E6055	MAS 200R Autosampler Piston O-Ring 29002195 Pack of 1	29002195
E3022	Nitrogen Column Flash 1112 0.5m 26007820	26007820
E4410	Replacement furnace Flash EA with 18mm and 25mm adaptors	35460018 35460021 35460022
F5112	Sealing Device 25209010	25209010
F5111	Sealing Device Solid Samples ANA1400 20500310	20500310
E2158	Silicone Grease 501-241 100gm	-
E2039	Snoop Leak Detection Fluid	-
F5104	Spatula (Left Handed) 20500600	20500600

F5103	Spatula 20500600	20500600
E6050	Top O-Ring For Alloy Reactor Flash EA1112 29020682 pack of 5	29020682
F4523	Tubing Nuts 2mm 35040403 pack of 10	1121370 35040403
E2040	Water/Carbon Dioxide Scrubber 28113060	28113060
E2047	Water/Carbon Dioxide Scrubber 28113106	28113106
E2059	Water/Carbon Dioxide Scrubber Glass Only Part of E2047	-
E2340	Water/CO2-trap for NA2000 Plastic part only (Complete trap with metal caps see E2040)	-

Prepacked reactors and consumables kits

Cat #	Description	OEM #
X1097	Consumable Kit 1000 Analyses N Org EA112 190 02413	19002413
X3211	Kit 1000 analyses NC Flash1112 1142960	1142960
X1150	Prepacked Oxidation High Performance Alloy Reactor 46802008	46802008

Quartz, Glassware and Furnace tubes

Cat #	Description	OEM #
C1100	Condenser Flash 1112	-
F5150	High Performance Alloy Reactor Flash EA1112 46820075	46820075
F5153	High Performance Inconel Reactor Crucible Flash EA1112 25205011	25205011
C1095	Quartz Ash Crucible Flash EA1112 N/Protein	-
C1092	Quartz Reactor(version of HPAR46820075) Not suitable for samples with high sodium content - Do not exceed 200kpa	-
C1087	Reduction Tube Quartz Flash Protein & Soils	-

Reference Materials

Cat #	Description	OEM #
B2140	Pasta Standard 33840024 10gm	33840024
B2158	Sorghum Flour Standard OAS 30gm See Certificate 114856 Expiry 16-Feb-12	-
B2142	Tryptophan OAS See Certificate For Values 100gm	-
B2130	Urea OAS See Certificate 2874 Expiry 02 jun 14 10gm	33840002
B2119	Urea OAS See Certificate For Values 25gm	-
B2156	Wheat Flour Standard OAS 30g See Certificate 114857 Expiry 16-Feb-12	-

Sample encapsulation

Cat #	Description	OEM #
F5109	Sealing Device for Liquids 20503002	20503002
D1078	Tin Discs 35mm pack of 100	25208018
D1066	Tin Foil Discs Standard Weight 30mm diameter pack of 100	25208015
D1059	Tin Foil Discs Standard Weight 34mm diameter pack of 100	-
D5055	Tin Weighing Pans Pressed Standard Weight 10 x 10mm pack of 100	25208000

Product Sheet ~ Thermo NA1112Protein

Placing your order: We are happy to accept your order by email, telephone, fax or post. We welcome all orders, large or small, with no minimum order requirement. Visa and Mastercard payments are accepted or we can open an account for you – please enquire. Most orders are despatched from stock within 24 hours.

Quality: All products are manufactured in our ISO9001:2000 approved facility and are covered by our 'Total Satisfaction Guarantee'.

Can't see what you need? We offer thousands of products for all elemental analysers. If you can't see what you need here the chances are we will have the product available. Just contact our friendly sales team who will track down what you need.

OEM Part numbers: The use of OEM part numbers is intended for convenience only and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and suitable for use in the stated instrument. All trade marks acknowledged.
